

POST-OPERATIVE  
GARMENTS AND BRAS



REVÉE®

ITALIAN COMPANY  
BASED IN TURIN  
SINCE 2012

## INDEX

---

04	POST-SURGICAL BRAS
20	BREAST AUGMENTATION BELT
24	LIPOSUCTION GARMENTS
38	BBL GARMENTS
44	TUMMY TUCK GARMENTS
54	REVÉE SLEEVES
60	FACE GARMENTS
64	REVÉE PADS
68	ANTITHROMBOSIS STOCKINGS
72	MEN'S GARMENTS
78	REVÉE CARDIO
86	REVÉE SCAR GEL
90	REVÉE SUNCARE 50+
92	REVÉE OIL
94	REVÉE SUPPORT
96	REVÉE NEWS
98	REFERENCES



REVÉE®



Company with **certified**  
**quality management**  
**system** TÜV Italia Srl  
UNI EN ISO 9001:2015  
**No. 50 100 16173 Rev.001**




REVÉE<sup>®</sup>  
your **body** our **care**

— [reeve.it](http://reeve.it) —



# POST-SURGICAL BRA

Our **Revée Bra** Collection offers a full solution for **post-operative recovery following breast augmentation, breast reduction, breast reconstruction**. It combines an exclusive design, engineered with the support of world reknown plastic surgeons with the highest-standards manufacturing process. Thanks to the unique characteristics of the **Lycra®** fabric (Latex free) it ensures optimal compression while providing maximum comfort to facillitate healing and repair of the tissues and infra-mammary structures (1). Patients who have undergone **breast augmentation surgery have more success and satisfaction with their procedure when they wear a quality compression post-surgical bra.**



*Learn more*

## AESTHETIC BREAST SURGERY

---


Revée Bra Laetitia  
REV.0301



## RECONSTRUCTIVE SURGERY

---


Revée Bra Julie  
REV.0303



## SECOND STAGE

---

Work out with Revée Bra Elodie  
REV.0302



# WE ARE REVEÉ.

Italian company specialised  
in the production of post-operative  
garments for Plastic, Reconstructive  
and Aesthetic Surgery.







*Learn more*

**Application fields:**  
breast augmentation,  
mastopexy, sternotomy,  
breast reduction, breast  
reconstruction.



Fully  
openable



Designed  
by surgeons



Ultra-flat  
seam



Expandable  
cups



No  
underwires



Breathable  
and quick  
to dry



3 front  
hooks



4 adjustable  
hooks



# LAETITIA BRA

REVÉE®

REV.0301



- 1 4 adjustable hooks on the shoulder straps.
- 2 Air Feel fabric on the back.
- 3 3 adjustable hooks on the front.
- 4 Extra breast support in Lycra®.
- 5 Soft jaquard elastic bend by Revée® of 4.5 cm.
- 6 External seam to avoid irritating the nipple and areola.
- 7 Bra cups self-mold and perfectly adapt themselves in case of breast asymmetry.


Product code: REV.0301

Size: XS - S - M - L - XL - XXL - 3XL - 4XL




Colours: 

Composition: 74% polyamide, 26% elastan

Washing: 



SIZE TABLE ON PAGE 16







# JULIE BRA

REVÈE®

**Application fields:**

breast augmentation, reduction, breast reconstruction and implants substitution, mastopexy.



Fully openable



Designed by surgeons



Dry fast



Expandable cups



No underwires



Ultra-flat seam



3 front hooks



4 adjustable hooks

# JULIE BRA

REVÉE®

REV.0303



- 1 3 adjustable hooks on the front.
- 2 Hook straps adjustable in 4 positions.
- 3 Additional X-shape support in Lycra® on the back.
- 4 Soft jaquard elastic bend by Revée® of 4.5 cm.
- 5 Bra cups self-mold and perfectly adapt themselves in case of breast asymmetry.


Product code: REV.0303

Size: XS - S - M - L - XL - XXL - 3XL - 4XL



Colours: 

Composition: 89% polyamide, 11% elasthan

Washing: 



SIZE TABLE ON PAGE 16







**ELODIE BRA**

REVÉE®

15



*Learn more*



Fully openable



Designed by surgeons



Dry fast



Expandable cups



No underwires



Ultra-flat seam



Front zip



4 adjustable hooks



# ELODIE BRA

REVÉE®

REV.0302



- 1 Closure: zipper and hooks.
- 2 Hook straps adjustable in 4 positions.
- 3 Additional X-shape support in Lycra® on the back.
- 4 Soft jaquard elastic bend by Revée® of 4.5 cm.
- 5 Bra cups self-mold and perfectly adapt themselves in case of breast asymmetry.


Product code: REV.0302

Size: XS - S - M - L - X - XXL - 3XL - 4XL




Colours: 


Composition: 89% polyamide, 11% elasthan

Washing: 



SIZE TABLE ON PAGE 16





LAETITIA  
REV.0301

ELODIE  
REV.0302

JULIE  
REV.0303

SIZE	UNDER BREAST	UNDER BREAST	UNDER BREAST
XS	66 - 70	67 - 72	67 - 72
S	71 - 75	73 - 77	73 - 77
M	76 - 80	78 - 82	78 - 82
L	81 - 85	83 - 87	83 - 87
XL	86 - 90	88 - 92	88 - 92
XXL	91 - 95	93 - 97	93 - 97
3XL	96 - 100	98 - 102	98 - 102
4XL	≥ 101	≥ 103	≥ 103

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

### ⚠ Warnings:

- In order to **guarantee its efficacy**, tolerability and proper use, **it should be worn with extreme care**.
- In case of concern on the application procedure, contact the shop staff or the treating physician.
- Choose the correct size that your specialist shop or treating physician recommends.
- The **fabrics** used are **hypoallergenic, non-cytotoxic, and latex-free**. It is however advisable to check the

components listed on the label to rule out any allergy to the fabrics used.

- In case a **surgical wound** is present, do not put into direct contact until complete healing.

■ **Close all the velcro and hook fasteners** before washing in order not to damage other items and to protect the shaping feature of the product.

**Tip:** For a better post-surgical recovery, acquiring two bras of the right size is recommended.

# REVÉE®

your **body** our **care**



# REVÉE BREAST AUGMENTATION BELT

## STANDARD



REV.0201

Revée Belt is a comfortable elastic band with an adjustable closure velcro to stabilize and avoid prosthesis mis-positioning. It is recommended after breast augmentation surgery in combination with Revée Bra.

**Application fields:** recommended after breast augmentation surgery or implants substitution.

## SHAPE



REV.0202

REV.0201



REV.0202





**REVÉE BELT**

REVÉE®



Enveloping  
belt



Designed  
by surgeons



Latex  
free



Adjustable  
velcro closure

## REV.0201



**BEST SELLER**

- 1 Adjustable velcro closure
- 2 Latex free
- 3 8-cm-high hook-and-loop containment elastic by Revée®

Product code: REV.0201

Size: One-size-fits-all

Colours:

Composition: 86% polyester, 14% rubber

Washing:

## REV.0202



- 1 Adjustable velcro closure
- 2 Latex free
- 3 Contoured under the armpits
- 4 Contoured under the armpits

Product code: REV.0202

Size: S - M - L


Colours:

Composition: 87% polyamide, 13% elastan

Washing:

SIZE	THORAX
S	70 - 85
M	86 - 105
L	> 105


\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS



## LIPOSUCTION GARMENTS


Revée® Lipo is a collection of multifunctional garments that meet all the specific requirements of an effective compression treatment, after liposuction interventions in one or more areas of the body: abdominoplasty, thighplasty, torsoplasty, body reshaping, lipofilling, and buttock lift.

Thanks to the **Lycra® fabric**, these garments are able to **self-mold** ensuring the best fit and **promoting the healing process**.



REV.0101 REV.0102 REV.0103





REV.0104   REV.0105   REV.0106

**REVÉE LIPO**

REVÉE

# REVÉE LIPO

REV.0101



- 1** Triple front hook closure
- 2** Soft jaquard elastic bend by Revée® of 4.5 cm.
- 3** Abdominal support
- 4** Wide hygienic opening


Product code: REV.0101 - 0102 - 0103

Size: XS - S - M - L - XL - XXL

Colours:

Composition: 74% polyamide,  
26% elasthan

Washing:



## REV.0102

## REV.0103

**BEST SELLER**



Abdominal support



Hygienic opening



3 front hooks



Designed by surgeons



Best seller

SIZE	WAIST	HIPS	THIGH	MID THIGH
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

# REVÉE LIPO

REV.0104



- 1 Triple front hook closure
- 2 Adjustable and removable straps
- 3 Abdominal support
- 4 Wide hygienic opening

Product code: REV.0104 - 0105 - 0106

Size: XS - S - M - L - XL - XXL

Colours:

Composition: 74% polyamide,  
26% elastan

Washing:



## REV.0105



**BEST SELLER**

## REV.0106



Abdominal support



Hygienic opening



3 front hooks



Designed by surgeons



Best seller

SIZE	WAIST	HIPS	THIGH	MID THIGH
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS




REVÉE®




# LIPOSUCTION GARMENTS

Revée® Lipo is a collection of multifunctional garments that meet all the specific requirements of an effective compression treatment, after liposuction interventions in one or more areas of the body: abdominoplasty, thighplasty, torsoplasty, body reshaping, lipofilling, and buttock lift. Thanks to the **Lycra® fabric**, these garments are able to **self-mold** ensuring the best fit and **promoting the healing process**.



Special liposuction support  
in the upper back



REV.0404 REV.0405 REV.0406

# REVÉE LIPO

REVÉE®

# REVÉE LIPO

REV.0404



- 1 Triple front hook closure
- 2 Adjustable and removable straps
- 3 Wide hygienic opening
- 4 Abdominal support
- 5 Additional back support in Lycra®



Product code: REV.0404 - 0405 - 0406

Size: XS - S - M - L - XL - XXL - 3XL

Colours:  

Composition: 74% polyamide,  
26% elastan

Washing:     





## REV.0405



## REV.0406



Abdominal support



Hygienic opening



3 front hooks



Designed by surgeons



Best seller

SIZE	WAIST	HIPS	THIGH	MID THIGH
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

# REVÉE LIPO

REV.0109

## Application fields:

Legs liposuction,  
postoperative care  
after a lipoedema



- 1 Triple front hook closure
- 2 Hygienic opening
- 3 Back seam
- 4 Foot extension

Product code: REV.0109


Size: XS - S - M - L - XL - XXL - 3XL - 4 XL

Height: Standard / Long


Colours: XXXX

Composition: 74% polyamide,  
26% elastan

Washing:



# REV.0109



Abdominal support



Hygienic opening



3 front hooks



Designed by surgeons




Best seller

SIZE	WAIST	HIPS	THIGH	MID THIGH
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66
4XL	117 - 129	132 - 148	89 - 99	62 - 70

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

# REVÉE LIPO



## Application fields:

Thigh lift, second post-liposuction phase.

- 1 Back seams to avoid friction in the scar area
- 2 Soft jaquard elastic bend by Revée® of 4.5 cm
- 3 Hook-free
- 4 Contoured garment

Product code: REV.0112

Size: XS - S - M - L - XL - XXL


Colours:

Composition: 74% polyamide, 26% elasthan


Washing:     

# REV.0112

## FRONT



## BACK



Designed  
by surgeons



Best  
seller




Latex  
free

SIZE	WAIST	HIPS	THIGH	MID THIGH
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS

# BBL GARMENTS

## STANDARD HOLES WITHOUT NET



REV.0107 REV.0108




*Learn more*

Revée® Lipo is a collection of multifunctional garments that meet all the specific requirements of an effective compression treatment, after liposuction interventions in one or more areas of the body: abdominoplasty, thighplasty, torsoplasty, body reshaping, lipofilling, and buttock lift.

Thanks to the **Lycra® fabric**, these garments are able to **self-mold** ensuring the best fit and **promoting the healing process**.



## BIGGER HOLES (BUTTOCKS&HIPS) WITH NET



REV.0407 REV.0408

# REVÉE BBL

REVÉE®

# REVÉE BBL



- 1 Adjustable and removable straps
- 2 Abdominal support
- 3 Closure: 3 positions of hooks
- 4 Opening on the buttocks
- 5 Wide hygienic opening


Product code: REV.0107 - 0108

Size: XS - S - M - L - XL - XXL - 3XL

Colours:

Composition: 74% polyamide,  
26% elastan


Washing:     






# STANDARD HOLES WITHOUT NET

REV.0107



REV.0108



Abdominal  
compression



3 front  
hooks



Designed  
by surgeons



Hygienic  
opening



Latex  
free



Best  
seller

SIZE	WAIST	HIPS	THIGH	MID THIGH
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS

# REVÉE BBL

(buttocks and hips)



- 1 Adjustable and removable straps
- 2 Abdominal support
- 3 Triple front hook closure
- 4 Bigger holes in the back (buttocks and hips)
- 5 Wide hygienic opening
- 6 Additional back support in Lycra®

Product code: REV.0407 - 0408

Size: XS - S - M - L - XL - XXL - 3XL


Colours:

Composition: 74% polyamide, 26% elasthan


Washing:     

# BIGGER HOLES WITHOUT NET

REV.0407



REV.0408



Abdominal  
compression



3 front  
hooks



Designed  
by surgeons



Hygienic  
opening



Latex  
free



Best  
seller

SIZE	WAIST	HIPS	THIGH	MID THIGH
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS

# TUMMY TUCK GARMENTS



BEST SELLER



REV.1351

After undergoing **abdominoplasty** surgery, it is essential to have **support in the post-operative** period. In collaboration with our **team of experts** we have created a line of compression garments specially developed to **improve your recovery** after abdominoplasty surgery and emphasise the postoperative results.

**Application fields:** post-operative compression therapy following an operation of: hip liposuction abdominoplasty, torsoplasty and lipofilling with harvesting in the abdominal region, caesarean section.



REV.1352

# REVÉE®

your **body** our **care**

REVÉE BODY

REVÉE®

45



Abdominal  
compression



Side  
zip



Designed  
by surgeons




Best  
seller



*Learn more*

# REVÉE BODY

REV.1351



- 1 Adjustable and removable straps
- 2 Abdominal support
- 3 Side closure

Product code: REV.1351

Size: XXS - XS - S - M - L - XL - XXL - 3XL

Colours:

Composition: 74% polyamide,  
26% elastan

Washing: 

SIZE	WAIST	HIPS
XXS	39 - 51	70 - 86
XS	49 - 61	78 - 94
S	59 - 71	86 - 102
M	69 - 81	94 - 110
L	79 - 91	102 - 118
XL	89 - 101	110 - 126
XXL	99 - 111	118 - 134
3XL	109 - 121	126 - 142

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS



← WAIST

← HIPS



Abdominal  
compression




3 front  
hooks



Designed  
by surgeons



Latex  
free



*Learn more*



# REVÉE BODY

## REV.1352



- 1 Adjustable and removable straps
- 2 Abdominal support
- 3 Front closure
- 4 Additional back support in Lycra®

Product code: REV.1352

Size: XS - S - M - L - XL - XXL - 3XL

Colours:

Composition: 74% polyamide,  
26% elastan

Washing: 

SIZE	WAIST	HIPS
XS	39 - 51	70 - 86
S	49 - 61	78 - 94
M	59 - 71	86 - 102
L	69 - 81	94 - 110
XL	79 - 91	102 - 118
XXL	89 - 101	110 - 126
3XL	99 - 111	118 - 134

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS

## REV.1353

**BEST SELLER**



- 1 Abdominal support
- 2 Closure: 3 positions of hooks on both sides
- 3 Soft jacquard elastic bend by Revée®

Product code: REV.1353

Size: XS - S - M - L - XL - XXL - 3XL - 4XL

Colours:

Composition: 74% polyamide,  
26% elastan

Washing:

SIZE	HIPS
XS	74 - 80
S	81 - 87
M	88 - 94
L	95 - 101
XL	102 - 108
XXL	109 - 115
3XL	116 - 122
4XL	123 - 129

**Application fields:** post-surgical compression therapy after surgeries such as: caesarean section, hysterectomy, mini-abdominoplasty, lipofilling with extraction from the lower abdominal region, and lower-abdominal surgeries in general.

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS



Abdominal compression



3 front hooks



Designed by surgeons



Ideal for post-partum



Adjustable on hips



*Learn more*



**Application fields:**

abdominoplasty,  
liposuction of the  
abdomen, hips,  
abdominal surgery,  
laparocoele. (2)\*



*Learn more*



Abdominal  
compression



Diastasis  
recti



Adjustable  
velcro closure



Best  
seller

# REVÉE BINDER

REVÉE®

REV.0203



cm26




**BEST SELLER**

REV.0204



cm32



- 1 Breathable fabric
- 2 Contoured
- 3 Easy to put on
- 4 Abdominal support
- 5 Front cotton  reinforcement

Product code: REV.0203 - 0204

Size: S - M - L - XL - XXL - 3XL

Colours:  





Height: 26 cm - 32 cm

Composition: 87% polyamide,  
13% elastan

Washing:     

SIZE	HIPS
S	70 - 82
M	82 - 94
L	94 - 105
XL	105 - 117
XXL	117 - 130
3XL	130 - 150

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS





Adjustable with  
hooks in the  
forearm area



Best  
seller




Designed  
by surgeons



**REVÉE SLEEVES**

REVÉE®

55



2 rows  
of hooks



Adjustable  
velcro closure



Best  
seller



Designed  
by surgeons




*Learn more*



# REVÉE SLEEVES

REV.200

**BEST SELLER**



- 1 Adjustable velcro closure
- 2 Closure: 2 positions of hooks

Product code: REV.200

Size: XS - S - M - L - XL - XXL

Colours:  

Composition: 74% polyamide,  
26% elastan

Washing:     

SIZE	ARM
XS	24 - 25.5
S	25.5 - 27
M	27 - 28.5
L	28.5 - 30
XL	30 - 34
XXL	≥ 35

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS

# REVÉE SLEEVES

REVÉE®

## REV.201



- 1 Upper back compression
- 2 4 rows of hooks for upper adjustment
- 3 Triple front hooks for abdominal adjustment

Product code: REV.201

Size: XS - S - M - L - XL - XXL

Colours:

Composition: 74% polyamide,  
26% elastan

Washing:

SIZE	ARM	UNDER BREAST
XS	24 - 26	64 - 69
S	26 - 27	69 - 74
M	27 - 29	74 - 79
L	29 - 30	79 - 84
XL	30 - 33	84 - 91
XXL	33 - 36	91 - 102

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

**Application fields:** post-surgical compression therapy after brachioplasty, arm lifting surgery, armpit surgery and back liposuction.



Adjustable hooks



Best seller



Designed by surgeons





Head-neck  
Revée  
garments



Adjustable  
velcro  
closure



Unisex

## FACE GARMENTS

Face compression garments provide targeted compression for chin liposuction or other facial surgical procedures. They are ideal for face and neck compression after neck lifts or chin liposuction. If applied correctly, these garments will increase blood circulation around the surgical area where

the surgery has been performed, as well as reduce swelling, eliminate any harmful bodily fluids, and speed up the healing process. In case of a sudden movement (e.g. sneeze or cough), these facial garments will minimize the discomfort.

# REVÉE MASKS

TOP OF HEAD-CHIN



## REV.5001



Product code: REV.5001

Size: S - M - L

Colours:  

Composition: 74% polyamide,  
26% elasthan

Washing:     

**Application fields:** lypoplasty  
of the chin, face lift.

SIZE	TOP OF HEAD-CHIN
S	53 - 57
M	58 - 62
L	63 - 67

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS

## REV.5002

**BEST SELLER**



Product code: REV.5002

Size: S - M - L - XL


Colours:  

Composition: 74% polyamide,  
26% elasthan

Washing:     

**Application fields:** forehead and face lift, lypoplasty of the chin, otoplasty, malar augmentation.

## REV.5003



Product code: REV.5003

Size: One-size-fits-all

Colours:  

Composition: 86% polyester,  
14% rubber

Washing:     

**Application fields:** lypoplasty of the chin, face lift, otoplasty.

SIZE	TOP OF HEAD-CHIN
S	53 - 57
M	58 - 62
L	63 - 67
XL	68 - 72



Unisex



Abdominal  
compression



Best  
seller



Cotton



Designed  
by  
surgeons



# REVÉE PAD

Revée Pads are worn under the garments to provide **greater compression on the treated area** following reconstructive, aesthetic or reconstructive surgery. Product available in **three different formats** with different intended uses. The pads provide a smooth and soft surface in contact with the skin, **improving comfort** after liposuction or after a body contouring surgery, **help**

**minimise the risk of haematomas.**

They also prevent the formation of bruises and skin folds, **optimising lymphatic drainage** by providing adequate compression and they help **accelerating the healing** by preventing water retention and reducing the accumulation of fluid in stubborn areas. Also indicated for patients suffering from lipoedema.



# REVÉE PAD

## REVÉE PAD01



UNISEX

Product code: REVÉE PAD01

Size: S - M / L - XL / XXL - 3XL

Colours:


Composition: 44% cotton,  
56% polyurethane

Washing:     

SIZE	WAIST
S - M	63 - 76
L - XL	77 - 94
XXL - 3XL	95 - 112

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS

## REVÉE PAD02



Product code: REVÉE PAD02

Size: One-size-fits-all

Colours:

Composition: 44% cotton,  
56% polyurethane

Washing:

## REVÉE PAD03



Product code: REVÉE PAD03

Size: One-size-fits-all

Colours:

Composition: 44% cotton,  
56% polyurethane

Washing:



2 pieces REVÉE PAD02  
in each pack



# Thrombosafe®

Graduated compression (more intense in the ankle, progressively decreasing towards the calf and thigh), which restores the blood circulation speed to normal levels in leg veins.

Thrombosafe® stockings stand out for their premium quality and reliability: gimped latex-free elastomer yarn covered with a protective sheath, for **excellent skin tolerability**. Elasticated edges with soft seams. The **inspection hole** allows the doctor to conveniently check the foot. Anatomic, seamless, enveloping and soft heel. Elastic non-slip band applied in the extremity preventing the stockings from going down and keeping them in the correct position. **LATEX-FREE** fabric, resistant to the use of skin creams, ointments, and sweat.

THROMBOSAFE®  
STOCKINGS



*Learn more*



Unisex




Latex free

REVEE®

# Thrombosafe®


REV.1609



Product code: REV.1609  
Size: S - M - L - XL (3 lengths)  
Colours:   
Composition: 80% Nylon  
20% Lycra  
Washing:

## How to choose your size

To choose the right size for the patient, measure the limb circumference (CB, CD, CG) to determine the stockings' width (SMALL, MEDIUM, LARGE, EXTRA LARGE). Then measure the length of the limb (iD or iG) to determine the stockings' length (short, standard, long). The colour facilitates the choice of the size and length. Twelve sizes available for each type of stockings, to ensure a correct fit for the patient's body.



S	Circumference cm			Length cm	
	cB	cD	cG	iD	iG
Short	19 - 22	28 - 34	45 - 55	34 - 39	65 - 75
Standard	19 - 22	28 - 34	45 - 55	39 - 44	75 - 85
Long	19 - 22	28 - 34	45 - 55	44 - 49	85 - 95

M	Circumference cm			Length cm	
	cB	cD	cG	iD	iG
Short	22 - 25	34 - 40	50 - 60	34 - 39	65 - 75
Standard	22 - 25	34 - 40	50 - 60	39 - 44	75 - 85
Long	22 - 25	34 - 40	50 - 60	44 - 49	85 - 95

L	Circumference cm			Length cm	
	cB	cD	cG	iD	iG
Short	25 - 28	40 - 46	55 - 65	34 - 39	65 - 75
Standard	25 - 28	40 - 46	55 - 65	39 - 44	75 - 85
Long	25 - 28	40 - 46	55 - 65	44 - 49	85 - 95

XL	Circumference cm			Length cm	
	cB	cD	cG	iD	iG
Short	28 - 32	46 - 53	60 - 75	34 - 39	65 - 75
Standard	28 - 32	46 - 53	60 - 75	39 - 44	75 - 85
Long	28 - 32	46 - 53	60 - 75	44 - 49	85 - 95

# GYNECOMASTIA GARMENTS GARMENTS

Our **Revée men's Collection** provides a complete solution for post-operative recovery following gynecomastia, lipoplasty of the back, midriff and hips, liposuction and abdominal surgery. It combines an exclusive design, engineered with the support of world reknown plastic surgeons with the highest-

standard manufacturing process. Thanks to the unique characteristics of the **Lycra® fabric** (Latex free) it ensures optimal compression while providing maximum comfort to facillitate healing and repair of the tissues and infra-mammary structures.



REV.030



REV.032





REV.031

REVÉE®

73

# REVÉE CHEST

REV.030






- 1 Two rows of front hooks
- 2 Pectoral support
- 3 Adjustable velcro closure

Product code: REV.030

Size: XXS - XS - S - M - L - XL

Colours:

Composition: 74% polyamide,  
26% elastan

Washing:     

SIZE	THORAX
XXS	72 - 80
XS	80 - 88
S	89 - 97
M	98 - 106
L	107 - 115
XL	116 - 124

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS





**Application fields:**

post-surgical compression therapy after the following surgeries: gynecomastia, surgical repair of the pectoral region, thoracic surgery or back liposuction.



*Learn more*



Pectoral compression velcro closure



Adjustable



3 front hooks



Best seller

# REVÉE CHEST

REV.031



- 1 Two rows of front hooks
- 2 Abdominal support
- 3 Pectoral support
- 4 Additional X-shape support
- 5 Velcro closure

Product code: REV.031

Size: XXS - XS - S - M - L - XL



Colours:

Composition: 74% polyamide, 26% elastan

Washing:     

SIZE	THORAX
XXS	72 - 80
XS	80 - 88
S	89 - 97
M	98 - 106
L	107 - 115
XL	116 - 124

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS



**Application fields:**

gynecomastia surgery,  
surgical repair of the male  
pectoral region, torsoplasty,  
flank liposuction, and  
abdominoplasty.



Pectoral  
compression



Adjustable  
velcro closure



3 front  
hooks

# REVÉE CARDIO

## REV.032



Revée Cardio can aid recovery after open heart surgery or other chest surgeries. It reduces pain, decreases the risk of complications, minimizes friction between the two halves of the sternum, while allowing to breathe more easily.

- 1 4 adjustable velcro closures
- 2 Sternal support
- 3 Contoured binder
- 4 Adjustable and removable straps



Product code: REV.032

Size: S - M - L - XL - XXL

Colours:           

Composition: 74% polyamide,  
26% elastan

Washing:

SIZE	THORAX
S	78 - 88
M	89 - 98
L	98 - 108
XL	109 - 118
XXL	119 - 128

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS

**Application fields:**

The garment has special handles that ensure maximum support when mostly needed. Coughing, deep breathing and returning to normal activities can be very painful after surgery. (3)\*



*Learn more*



Pectoral compression velcro closure



Adjustable velcro closure





Latex free



Ideal after sternotomy procedures

# REVÉE LIPO FOR MAN

REV.0104



SIZE	WAIST	HIPS	THIGH	MID THIGH
XS	47 - 59	76 - 92	40 - 50	34 - 42
S	57 - 69	84 - 100	47 - 57	38 - 46
M	67 - 79	92 - 108	54 - 64	42 - 50
L	77 - 89	100 - 116	61 - 71	46 - 54
XL	87 - 99	108 - 124	68 - 78	50 - 58
XXL	97 - 109	116 - 132	75 - 85	54 - 62
3XL	107 - 119	124 - 140	82 - 92	58 - 66

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS



Our **Revée men's Collection** provides a complete solution for post-operative recovery following gynecomastia, lipoplasty of the back, midriff and hips, liposuction and abdominal surgery. It combines an exclusive design, engineered with the support of world reknown plastic surgeons with the highest-standard manufacturing process. Thanks to the unique characteristics of the Lycra® fabric (Latex free) it ensures optimal compression while providing maximum comfort to facillitate healing and repair of the tissues and inframammary structures.



- 1 Closure: 3 positions of hooks
- 2 Adjustable and removable straps
- 3 Abdominal support

Product code: REV.0104

Size: XS - S - M - L - XL - XXL

Colours:

Composition: 74% polyamide,  
26% elasthan

Washing:

**Application fields:** Abdominoplasty, liposuction of the abdomen, hips, thighs and lipofilling.




Designed  
by  
surgeons



Abdominal  
compression

# REVÉE LIPO FOR MAN

REV.0504



SIZE	THORAX	WAIST	HIPS	THIGH
XS	87 - 91	75 - 79	82 - 85	47 - 50
S	92 - 96	80 - 86	86 - 90	51 - 55
M	97 - 102	87 - 93	91 - 97	56 - 60
L	103 - 109	94 - 101	98 - 104	61 - 65
XL	110 - 117	102 - 110	105 - 112	66 - 70
XXL	118 - 125	111 - 119	113 - 120	71 - 75

\*NOTES: MEASUREMENTS ARE INDICATIVE AND EXPRESSED IN CENTIMETERS

- 1 Closure: 3 positions of hooks
- 2 Abdominal support
- 3 Chest support
- 4 Shoulder closure with three rows of hooks

Product code: REV.0504

Size: XS - S - M - L - XL - XXL

Colours:

Composition: 74% polyamide,  
26% elastan

Washing:     

**Application fields:** gynecomastia surgery, reconstruction of the male pectoral region, trochanteric region (upper thighs), hip liposuction, abdominoplasty.



Pectoral support



Designed by surgeons



Abdominal compression



### Application fields:

Abdominoplasty,  
liposuction of the  
abdomen, hips, abdominal  
surgery, laparocoele.(2)\*



Abdominal  
compression



Adjustable  
velcro closure



Best  
seller

# REVÉE BINDER

REV.0203



REV.0204



- 1 Breathable fabric
- 2 Contoured
- 3 Easy to put on
- 4 Abdominal support
- 5 Front cotton reinforcement

Product code: REV.0203 - 0204

Size: S - M - L - XL - XXL - 3XL

Colours:

Height: 26 cm - 32 cm

Composition: 87% polyamide,  
13% elastan

Washing:

SIZE	HIPS
S	70 - 82
M	82 - 94
L	94 - 105
XL	105 - 117
XXL	117 - 130
3XL	130 - 150

\*NOTES: MEASUREMENTS ARE INDICATIVE  
AND EXPRESSED IN CENTIMETERS



Available in three **formats**:  
**10 g - 20 g** and **60 g**



Medical silicone  
Odourless  
Without parabens



*Learn more*



Silicone gels are the first treatment recommended by international guidelines to cure and prevent hypertrophic scars and keloids. (4)\*

# REVÉE®

## SCAR GEL



IN THE BEST PHARMACIES AND HEALTHCARE SHOPS

# REVÉE®

---

## SCAR GEL



It reduces redness and itching, while improving the aesthetical aspect. The scar becomes smoother and more elastic.



Odourless.



Quick-drying silicone protective film.

## TREATMENT FOR SCARS RECOMMENDED BY PLASTIC SURGEONS

It is indicated for the treatment and prevention of hypertrophic or keloid scars caused by a burn, a surgical procedure or trauma.

It is also used for immature and red scars in relief during the active healing phase.



### How to use:


- 1 Clean the area to be treated. Then dry by dabbing;
- 2 Apply a thin layer of Revée Scar Gel. Do not massage;
- 3 After Revée® Scar Gel is dried you can apply: makeup, cosmetics or sunscreen;
- 4 Apply twice a day, morning and evening for at least 2/3 months.

**Ingredients:** Cyclopentasiloxane, cyclohexasiloxane, dimethicone/vinyl dimethicone crosspolymer.



## THYROIDECTOMY

scar




## TRAUMA

scar




## MASTOPLASTY

scar



## NEVUS REMOVAL

scar



## C-SECTION and ABDOMINOPLASTY

scar





In the best pharmacies



Very water resistant  
**High protection**  
Antioxidant  
Non-perfumed cream  
Nourishing  
Nickel tested  
Gluten free

75ml

REVEE<sup>®</sup>  
**50+**  
sun/care



GLUTEN FREE

Sunscreen SPF 50+.  
**DERMATOLOGICALLY TESTED.**  
Highly resistant to water. High protection

[www.revee.it](http://www.revee.it)



8 051122 405305




UVB + UVA

# Protect Scars from the Sun.




*Learn more*



# DRY BASE OIL WITH

- ✓ Elasticising
- ✓ Emollient
- ✓ Soothing
- ✓ Hydrating
- ✓ Protective



+ SWEET ALMOND OIL



+ ROSE HIP



+ BIO-FERMENTED  
VEGETABLE OILS



+PURE VITAMINE E



It improves and reduces stretch marks.  
It counteracts skin ageing and dehydration



REVÉE®

tiO



NON COMEDOGENIC

**NO**  
FRAGRANCES  
PRESERVATIVES  
DYES

#REVÉESUPPORT



REVÉE®

Revée proudly supports research through the project

# SUTURE IN SPACE


#SUTUREINSPACE



Video YouTube



Learn more




**Melanoma: new approaches and implications for prevention**

Cutaneous melanoma represents a growing challenge for the medical community, with a significant increase in cases over the last decades....

MEDICAL DOCTOR AND SURGEONS



**Breast Units: a multidisciplinary approach against breast cancer**



**LESC: the innovative approach for reducing localised adiposity**



Now also on **WhatsApp!**



REVÉE<sup>®</sup>  
**NEWS**

discover our newspaper

[www.revee.news](http://www.revee.news)



## REFERENCES

- 1 - Chen et al. 1998 ***“The effect of a pressure garment on post-surgical drainage and seroma formation in breast cancer patients.”*** Singapore Med J. 1998 Sep;39(9):412-5.
- 2 - Rothman et al. 2014 ***“Abdominal binders may reduce pain and improve physical function after major abdominal surgery - a systematic review.”*** Dan Med J. Nov;61(11):A4941.
- 3 - Klement et al. 2010 ***“Supportive treatment using a compression garment vest of painful sternal instability following deep surgical wound infection: a case report.”*** J Med Case Rep. 2010 Aug 11;4:266. doi: 10.1186/1752-1947-4-266.
- 4 - Hsu et al. 2017 ***“Review of Silicone Gel Sheeting and Silicone Gel for the Prevention of Hypertrophic Scars and Keloids.”*** Wounds. 2017 May;29(5):154-158.



R E V É E<sup>®</sup>

[www.revee.it](http://www.revee.it)

REVEE®  
your **body** our **care**

REVEE srl Via Tetti Castagno, 6 10020 Andezeno (TO) Italy

DEALER