



### Bowl Stands, Single-Double: 24555/24575

- Made of AISI 304 stainless steel
- Demontable body and 8 L/L bowls
- Ergonomic handle for easy moving
- 4 wheel profile feet



#### Technical Specifications

| Models       | Order Codes |
|--------------|-------------|
| 24555 Single | 14-5276     |
| 24575 Double | 14-8432     |



### Medical Waste Container: 24500/24510

- Made of stainless steel to be used in operating and intervention rooms
- 20 L internal capacity (24500)
- 50 L internal capacity (24510)
- Lid and carriage handle
- Rubber coated metal bumper against impacts



#### Technical Specifications

| Models         | Order Codes |
|----------------|-------------|
| 24500 20 Liter | 14-8897     |
| 24510 50 Liter | 14-7223     |



### Foot Step: 24015/24020/24025

- Made of stainless steel
- Non-slip, wide surface for easy use



#### Technical Specifications

| Models            | Dimensions (W x L x H) | Order Codes |
|-------------------|------------------------|-------------|
| 24015 Single Step | 290x500x180 mm         | 14-7322     |
| 24025 Double Step | 600x500x390 mm         | 14-7654     |



### Laundry Collection Trolley: 24760/24790

- Made of stainless steel
- Detachable and washable fabric collection bags
- Four castors, two of them with brake



#### Technical Specifications

| Models       | Dimensions (W x L x H) | Order Codes |
|--------------|------------------------|-------------|
| 24760 Single | 500x550x1000 mm        | 74-2452     |
| 24790 Double | 500x1100x970 mm        | 14-6082     |